

medi tech
systems limited

Supporting every breath

Electrosurgery and theatre consumables




INTERACT WITH US

facebook: meditechsystemslimited

twitter: @MeditechLtd

QUALITY | VALUE | INTEGRITY



18,000sq ft warehouse space.



Injection moulding processes.



Cleanroom facilities.

Meditech Systems Limited supply superior, affordable, quality accessories and theatre consumables.

Meditech is a well established British based manufacturer with a strong commitment to product innovation. This innovation is borne out of our close customer contact and has proven to be the catalyst for our growth.

Our comprehensive range covers all major makes of electrosurgery generators worldwide and is designed to meet the high quality and safety standards of each market. Our ability to provide the right product by design for specific clinical needs is nurtured by close liaison with our clinical partners. This approach provides benefits to both patient and clinician.

Our commitment is backed by increased capacity in our 18,000 sq ft UK based facility offering the following facilities: cleanroom for processing, in-house tooling design, injection moulding both automatic and insert moulding, extrusion and corrugation lines producing a variety of tubing sizes plus a nylon coating plant for insulation of our instruments.

Our aim at Meditech Systems Limited is to grow through partnership with our customers providing the quality, value and service required in today's ever changing health care industry.



SAVE MONEY on your clinic consumable purchasing.

To help you change, we will reduce your existing supplier prices by at least 5%

Do you want effective stock control which saves you money. We can offer a quick cost-effective solution with our consignment stock plan.

How it works

1. Stock holding agreed with the clinic staff.
2. Trolley and initial stock delivered free of charge by Meditech. (Please note the trolley and initial stock remain property of Meditech)
3. Our representative will visit every 4-6 weeks (or as agreed) to count the consignment stock.
4. Stock used is replenished and then charged for. The clinic has only to pay for items used.

Benefits

- You only pay for stock you have used.
- Stock replenished in singles, no minimum order quantity or box quantities.
- 30 days credit to pay for items used.
- Out of date stock? This would be our responsibility.
- No delivery charges.

CALL OUR TEAM
+44(0)1747 821546

Product coding formats

31 24 130 N 09

PENCIL ELECTRODES

- Digits relate to the product e.g. pencil electrodes.
- Digits relate to shaft size (2.4mm)
- Digits relate to length (130mm)
- Digits relate to needle style.
e.g. S-spatula N-needle B-ball HD-hook dagger
- Digits relate to tip diameter (0.9mm)

10 40 80 48 11 43 80 75

CABLES

- Digits relate to the product.
- Digits relate to generator connection/pin size.
- Digits relate to cable length (metres)
- Digits relate to instrument connection/socket style.

7 7 10 180 UK

FORCEPS

- Digits relate to the product e.g. bipolar forceps.
- Digits relate to forceps style.
- Digits relate to tip size.
- Digits relate to overall length.
- Letters relate to cable connection.

30 24 130 L 15

LOOP/BALL ELECTRODES

- Digits relate to the product.
- Digits relate to shaft diameter (2.4mm)
- Digits relate to overall length (130mm)
- Digits relate to style frame.
e.g. L-standard LP-low profile SQ-square T-ribar B-ball SA-sail
- Digits relate to loop diameter (15mm)

INDEX

Electrosurgery

- 06 Hand control pencils
- 07 Pencils, handles and electrodes
- 08 Electrodes
- 10 Patient plates
- 11 Monopolar cables
- 14 Bipolar cables
- 17 Bipolar forceps
- 20 Monopolar scissors
- 21 ENT instruments

Gynaecology

- 22 Loop and ball electrodes
- 24 Speculums
- 25 Biopsy instruments
- 26 Disposable consumables

Theatre & clinic consumables

- 27 Theatre and clinic consumable
- 28 Safety scalpels
- 29 Podiatry
- 30 Skin biopsy punch and dermal curettes
- 31 Thudicums and suction products

Berchold, Erbe, Martin, Eschmann, Valleylab, Bovie, Downs, Siemens, Wolf, Olsen, Storz, Martin, EMC, Dolley and Lamsley are brand names of respective companies




FDA REGISTRATION NO.10048739 AUSTRALIA TGA LISTED

Safety information


All Medtech products are designed and manufactured to meet the requirements of relevant BS, ISO and EN standards. Compliance with ISO 9001: 2008 ISO 13485 and EU Directive 93/42/EEC Annex II, Annex V of the Medical Devices Directive allows CE marking of all of our products. This enables products to be traded throughout Europe and the world without restriction.



Hand control pencils



Tip cleaner for cleaning electrodes. Code: 317520



| Product | Length | Code | Usage | Box qty | Product | Code |
|---------------------------------------|----------|-----------|----------------|---------|---|---------|
| 'DURA-SWITCH' REUSABLE PENCILS | | | | | | |
| Hand control pencil | 3 metres | 3103024 | 50 | 5 | Eshmann | 3144440 |
| Hand control pencil | 5 metres | 3105024 | 50 | 5 | Erbe | 3144450 |
| Hand control pencil | 3 metres | 3113024 | 25 | 5 | Neomed | 3144441 |
| Hand control pencil | 5 metres | 3115024 | 25 | 5 | Martin | 3144442 |
| Hand control pencil | 3 metres | SU3113024 | Single use (2) | 50 | Adaptors allow connection of hand control pencils (3 pin plug) to listed generators. Codes below show 3 and 5 metre length. | |
| Hand control pencil | 5 metres | SU3115024 | Single use (2) | 50 | | |
| ADAPTORS | | | | | | |

Generator connection plugs

Twin pedals



Single yellow



Single blue



Single black



Single white



FOOT PEDALS

end view

side view

| | | |
|--|--|-------------------------------------|
| | | Valleylab Conmed/Bardarm Erbe |
| | | Eschmann |
| | | GU-solar |
| | | Downs Diadon 450P2 |
| | | Sorling |

| | | | | |
|----------|----------|----------|----------|----------|
| 10313020 | 1031301Y | 1031301B | 10313010 | 1031301W |
| 10323020 | 1032301Y | 1032301B | 10323010 | 1032301W |
| 10343020 | 1034301Y | 1034301B | 10343010 | 1034301W |
| 10363020 | 1036301Y | 1036301B | 10363010 | 1036301W |
| 10373020 | 1037301Y | 1037301B | 10373010 | 1037301W |
| 10333020 | 1031301Y | 1031301B | 10313010 | 1031301W |
| 10383020 | 1031301Y | 1031301B | 10313010 | 1031301W |

Pencils, handles and electrodes

| Product | Length | Code | SU code | Connections |
|--------------------------|----------|----------|------------|-------------|
| SLIM-LINE PENCILS | | | | |
| Slimline pencil | 3 metres | 31403024 | SU31403024 | |
| Slimline pencil | 5 metres | 31405024 | SU31405024 | |
| Slimline pencil | 3 metres | 31803024 | SU31803024 | |
| Slimline pencil | 5 metres | 31805024 | SU31805024 | |
| Slimline pencil | 3 metres | 31503024 | SU31503024 | |
| Slimline pencil | 5 metres | 31505024 | SU31505024 | |

Handles, needles and electrodes



Screw fit 9BA electrodes

| Product | Code | Box qty |
|--------------------------|----------|---------|
| 9BA HANDLES | | |
| Bristol handle 120mm 9BA | 31481209 | 1 |
| 150mm 9BA | 31481509 | 5 |
| 200mm 9BA | 31482009 | 5 |

9BA THREAD ELECTRODES




| | | |
|-----------------------|---------|---|
| Thread electrodes 9BA | 3195506 | 5 |
| Thread electrodes 9BA | 3195508 | 5 |
| Thread electrodes 9BA | 3195510 | 5 |
| Thread electrodes 9BA | 3195512 | 5 |

9BA THREAD NEEDLES



| | | |
|--------------------|----------|---|
| Thread needles 9BA | 3195507 | 5 |
| Thread needles 9BA | 3193207 | 5 |
| Thread needles 9BA | 3195007A | 5 |
| Thread needles 9BA | 3195010 | 5 |



Slimline pencils with pushfit electrode connection. 2.4mm and 4mm sockets

| Product | Length | Dia | Code | Box qty |
|---|--------|-------|------------|---------|
| TUNGSTEN NEEDLES (COLORADO STYLE) | | | | |
| Needle Angled, insulated-exposed tungsten tip | 150mm | 0.7mm | 3624150N7A | 5 |
| Needle insulated-exposed tungsten tip | 70mm | 0.7mm | 3624070N7 | 5 |
| Needle insulated-exposed tungsten tip | 100mm | 0.7mm | 3624100N7 | 5 |
| Needle insulated-exposed tungsten tip | 150mm | 0.7mm | 3624150N7 | 5 |

0.5mm dia also available. For Single use prefix code with SU (box quantity 10)

Electrodes



4mm diameter 2.4mm

Product

Codes

2.4mm

30

4mm

30

2.4mm reusable

BLADE ELECTRODES

Blade electrode 2.4/4mm shaft 70mm length

SU3124070S

SU3140070S

3124070S

Blade electrode 2.4/4mm shaft 100mm length

SU3124100S

SU3140100S

3124100S

Blade electrode 2.4/4mm shaft 130mm length

SU3124130S

SU3140130S

3124130S

Blade electrode 2.4/4mm shaft 150mm length

SU3124150S

SU3140150S

3124150S

Blade electrode 2.4/4mm shaft 200mm length

SU3124200S

SU3140200S

3124200S



Product

Codes

2.4mm

30

4mm

30

2.4mm reusable

BLADE ANGLED ELECTRODES

Blade angled electrode 2.4/4mm shaft 70mm length

SU3124070A

SU3140070A

3124070A

Blade angled electrode 2.4/4mm shaft 100mm length

SU3124100A

SU3140100A

3124100A

Blade angled electrode 2.4/4mm shaft 130mm length

SU3124130A

SU3140130A

3124130A

Blade angled electrode 2.4/4mm shaft 150mm length

SU3124150A

SU3140150A

3124150A

Blade angled electrode 2.4/4mm shaft 200mm length

SU3124200A

SU3140200A

3124200A



Product

Codes

2.4mm

30

4mm

30

2.4mm reusable

NEEDLE ELECTRODES

Needle electrode 2.4/4mm shaft 70mm length 0.9 tip

SU3124070N09

SU3140070N09

3124070N09

Needle electrode 2.4/4mm shaft 100mm length 0.9 tip

SU3124100N09

SU3140100N09

3124100N09

Needle electrode 2.4/4mm shaft 130mm length 0.9 tip

SU3124130N09

SU3140130N09

3124130N09

Needle electrode 2.4/4mm shaft 150mm length 0.9 tip

SU3124150N09

SU3140150N09

3124150N09

Needle electrode 2.4/4mm shaft 200mm length 0.9 tip

SU3124200N09

SU3140200N09

3124200N09



Product

Codes

2.4mm

30

4mm

30

2.4mm reusable

NEEDLE ANGLED ELECTRODES

Needle angled electrode 2.4/4mm shaft 70mm length 0.9 tip

SU3124070N09A

SU3140070N09A

3124070N09A

Needle angled electrode 2.4/4mm shaft 100mm length 0.9 tip

SU3124100N09A

SU3140100N09A

3124100N09A

Needle angled electrode 2.4/4mm shaft 130mm length 0.9 tip

SU3124130N09A

SU3140130N09A

3124130N09A

Needle angled electrode 2.4/4mm shaft 150mm length 0.9 tip

SU3124150N09A

SU3140150N09A

3124150N09A

Electrodes



Product

Codes

2.4mm 20 4mm 20 3-8mm ball ends

BALL ELECTRODES

Ball electrode 2.4/4mm shaft 70mm length

SU3124070B4

SU3124070B4

Ball electrode 2.4/4mm shaft 100mm length

SU3124100B4

SU3124100B4

Ball electrode 2.4/4mm shaft 130mm length

SU3124130B4

SU3124130B4

Ball electrode 2.4/4mm shaft 150mm length

SU3124150B4

SU3124150B4

Ball electrode 2.4/4mm shaft 200mm length

SU3124200B4

SU3124200B4



Available in 3-8mm balls

Suffix all codes with your ball size requirements, i.e. 3124070B3, 3124070B4, 3124070B5, 3124070B6 etc.



Product

Codes

2.4mm reusable single use 20 electrode style

LAPAROSCOPIC ELECTRODES

Laparoscopic Hook electrode

3124200H

SU3124200H

Laparoscopic Hook and Dagger electrode

3124200HD

SU3124200HD

Laparoscopic Ball electrode

3840203A

SU3840203A



Product

Codes

2.4mm reusable single use 20

ANTHROSCOPIC HOOK - DAGGER ELECTRODES

Hook/dagger electrode 2.4mm shaft 180mm length

3124180HD

SU3324180HD

Hook/dagger electrode 2.4mm shaft 140mm length

3124140HD

SU3324140HD

Hook/dagger electrode 2.4mm shaft 120mm length

3124120HD

SU3324120HD



Product

Codes

2.4mm reusable single use 20

ANTHROSCOPIC HOOK ELECTRODES

Hook electrode 2.4mm shaft 180mm length

3124180H

SU3324180H

Hook electrode 2.4mm shaft 140mm length

3124140H


SU3324140H

Hook electrode 2.4mm shaft 120mm length

3124120H

SU3324120H

Patient plates



Equipotential ring design



Product

Codes

PATIENT PLATES

| | |
|---|-----------|
| Single plate adult | 901001 |
| Single plate paediatric | 901002 |
| Single plate adult and cable-Valleylab | 90103001 |
| Single plate paediatric and cable-Valleylab | 90103001P |
| Patient plate with 3m cable plug VL no centre pin | 90103002 |
| Split adult plate square shape butterfly | 901101 |
| Standard split plate adult | 901102 |
| Oval plate standard universal | 901210 |
| Oval split plate universal | 901220 |
| Child split plate | 901222 |
| Paediatric split plate | 901223 |
| Patient plate cable VL style no pin 3m length 25mm clamp | 902203025 |
| Patient plate cable VL style no pin 5m length 25mm clamp | 902205025 |
| Patient plate cable VL plug 5M length 25mm clamp | 902215025 |
| Patient plate adult with cable plus 4mm covered plug | 9040S30AD |
| Patient plate cable fixed 4mm pins 3m 25mm clamp | 904403025 |
| Dispersive cable 3M with twin plugs and patient plate clamp | 9044L3025 |
| Patient plate adult with cable 4m plug | 90601 |
| Patient plate infant with cable 4m plug | 90602 |
| Patient plate cable 6.4mm jack plug 3m 25mm clamp | 906403025 |
| Patient plate cable 6.4mm jack plug 5m 25mm clamp | 906405025 |

Patient plate clamp connection

Generator


Cable length

Code



PATIENT PLATE CABLES FOR VALLEYLAB GENERATOR CONNECTION

| | | |
|------------------|----------|-----------|
| Valleylab | 3 metres | 902213025 |
| Valleylab | 5 metres | 902215025 |
| Valleylab no pin | 3 metres | 902203025 |
| Valleylab no pin | 5 metres | 902205025 |



PATIENT PLATE CABLES FOR 6.4MM JACK PLUG GENERATOR CONNECTION

| | | |
|-----------------|----------|-----------|
| 6.4mm jack plug | 3 metres | 906403025 |
| 6.4mm jack plug | 5 metres | 906405025 |



PATIENT PLATE CABLES FOR TWIN 4MM PIN GENERATOR CONNECTION

| | | |
|--------------|----------|-----------|
| Twin 4mm pin | 3 metres | 904403025 |
| Twin 4mm pin | 5 metres | 904415025 |



Monopolar cables



4mm pin connection can be used with the following generators:
Berchtold, Erbe, Martin and Eschmann

| Product | Length | Code reusable | Code single use | Ends |
|--|----------|---------------|-----------------|------|
| 4mm PIN CONNECTION | | | | |
| Standard cable with socket 4.8 diameter | 3 metres | 10403048 | SU10403048 | |
| Standard cable with socket 4.8 diameter | 5 metres | 10405048 | SU10405048 | |
| Standard cable with socket 4.0 diameter | 3 metres | 10403040 | SU10403040 | |
| Standard cable with socket 4.0 diameter | 5 metres | 10405040 | SU10405040 | |
| Standard cable with socket 3.0 diameter | 3 metres | 10403030 | SU10403030 | |
| Standard cable with socket 3.0 diameter | 5 metres | 10405030 | SU10405030 | |
| Snare cable with socket 2.8-3.0 diameter | 3 metres | 10403028 | SU10403028 | |
| Snare cable with socket 2.8-3.0 diameter | 5 metres | 10405028 | SU10405028 | |
| Wolf coaxil cable with socket Pin 1.95 | 3 metres | 104030195 | SU104030195 | |
| Wolf coaxil cable with socket Pin 1.95 | 5 metres | 104050195 | SU104050195 | |
| Olsen coaxil cable | 3 metres | 10403042 | SU10403042 | |
| Olsen coaxil cable | 5 metres | 10405042 | SU10405042 | |
| Coaxil resectorscope (Ad/Pd) universal | 3 metres | 1040301U | SU1040301U | |
| Coaxil resectorscope (Ad/Pd) universal | 5 metres | 1040501U | SU1040501U | |
| Coaxil cable with socket pin 4mm | 3 metres | 10403044 | SU10403044 | |
| Coaxil cable with socket pin 4mm | 5 metres | 10405044 | SU10405044 | |

MEDI COUNT

NO MORE PUNCHING OUT HOLES OR TIME CONSUMING RECORD CARDS:

- **COPYRIGHT SOFTWARE:** records product details to RF datachip.
- **FACTORY BASE UNIT:** product details are added to the chip via computer, including usage and count setting (product repair history is recorded)
- **HOSPITAL UNIT RECORDS:** once processed simply insert the chip and automatic count is activated.
- **MICRO SIZE AUTOCLAVEABLE RF CHIP:** can be attached to a wide range of products.

LCD display indicates LIFE USES REMAINING and END OF LIFE FLASHING WARNING.



Monopolar cables



8mm Bovie connection can be used with the following generators: Valleylab, Erbe and Bovie

| Product | Length | Code reusable 1 | Code single use 10 2 | Ends |
|---------|--------|---------------------------------|---|------|
|---------|--------|---------------------------------|---|------|

8mm BOVIE CONNECTION

| | | | | |
|--|----------|-----------|-------------|--|
| Standard cable with socket 4.8 diameter | 3 metres | 10803048 | SU10803048 | |
| Standard cable with socket 4.8 diameter | 5 metres | 10805048 | SU10805048 | |
| Standard cable with socket 4.0 diameter | 3 metres | 10803040 | SU10803040 | |
| Standard cable with socket 4.0 diameter | 5 metres | 10805040 | SU10805040 | |
| Standard cable with socket 3.0 diameter | 3 metres | 10803030 | SU10803030 | |
| Standard cable with socket 3.0 diameter | 5 metres | 10805030 | SU10805030 | |
| Snare cable with socket 2.8-3.0 diameter | 3 metres | 10803028 | SU10803028 | |
| Snare cable with socket 2.8-3.0 diameter | 5 metres | 10805028 | SU10805028 | |
| Wolf coaxil cable with socket pin 1.95 | 3 metres | 108030195 | SU108030195 | |
| Wolf coaxil cable with socket pin 1.95 | 5 metres | 108050195 | SU108050195 | |
| Olsen coaxil cable | 3 metres | 10803042 | SU10803042 | |
| Olsen coaxil cable | 5 metres | 10805042 | SU10805042 | |
| Coaxil resectorscope (Ad/Pd) universal | 3 metres | 1080301U | SU1080301U | |
| Coaxil resectorscope (Ad/Pd) universal | 5 metres | 1080501U | SU1080501U | |
| Coaxil cable with socket pin 4mm | 3 metres | 10803044 | SU10803044 | |
| Coaxil cable with socket pin 4mm | 5 metres | 10805044 | SU10805044 | |

Generator adaptors



8mm Bovie connection with 4mm socket. Code: 108040

Monopolar cables



5mm pin connection can be used with the following generators:
Berchtold, Erbe, Martin and Eshmann

| Product | Length | Code reusable | Code single use | Ends |
|---------|--------|------------------|--------------------|------|
|---------|--------|------------------|--------------------|------|

5mm PIN CONNECTION

| | | | |
|--|----------|-----------|-------------|
| Standard cable with socket 4.8 diameter | 3 metres | 10503048 | SU10503048 |
| Standard cable with socket 4.8 diameter | 5 metres | 10505048 | SU10505048 |
| Standard cable with socket 4.0 diameter | 3 metres | 10503040 | SU10503040 |
| Standard cable with socket 4.0 diameter | 5 metres | 10505040 | SU10505040 |
| Standard cable with socket 3.0 diameter | 3 metres | 10503030 | SU10503030 |
| Standard cable with socket 3.0 diameter | 5 metres | 10505030 | SU10505030 |
| Snare cable with socket 2.8-3.0 diameter | 3 metres | 10503028 | SU10503028 |
| Snare cable with socket 2.8-3.0 diameter | 5 metres | 10505028 | SU10505028 |
| Wolf coaxil cable with socket pin 1.95 | 3 metres | 105030195 | SU105030195 |
| Wolf coaxil cable with socket pin 1.95 | 5 metres | 105050195 | SU105050195 |
| Olsen coaxil cable | 3 metres | 10503042 | SU10503042 |
| Olsen coaxil cable | 5 metres | 10505042 | SU10505042 |
| Coaxil resectorscope (Ad/Pd) universal | 3 metres | 1050301U | SU1050301U |
| Coaxil resectorscope (Ad/Pd) universal | 5 metres | 1050501U | SU1050501U |
| Coaxil cable with socket pin 4mm | 3 metres | 10503044 | SU10503044 |
| Coaxil cable with socket pin 4mm | 5 metres | 10505044 | SU10505044 |



CLICK COUNTER



THE CLICK COUNTER OFFERS THE FOLLOWING:

- **ADVANCE MARKER:** one click each time product is processed, it cannot be reversed.
- **ACCURATE RECORD:** easily visible scale shows accurate processing record.
- **FLEXIBLE SYSTEM:** easily attached to a wide range of devices.

Click counter
code: 17777

100

WE CAN ALSO OFFER TRADITIONAL PUNCH TRACKING CARDS.

Bipolar cables twin 4mm pin



Twin 4mm pins compatible with most generators

| Product | Length | Code reusable | Code single use | Ends |
|---------|--------|------------------|--------------------|------|
|---------|--------|------------------|--------------------|------|

TWIN 4mm PINS

| Product | Length | Code reusable | Code single use | Ends |
|----------------|----------|------------------|--------------------|------|
| Euro | 3 metres | 11443075 | SU11443075 | |
| Euro | 5 metres | 11445075 | SU11445075 | |
| Storz | 3 metres | 11443041 | SU11443041 | |
| Storz | 5 metres | 11445041 | SU11445041 | |
| Wolf | 3 metres | 11443042 | SU11443042 | |
| Wolf | 5 metres | 11445042 | SU11445042 | |
| US | 3 metres | 11443076 | SU11443076 | |
| US | 5 metres | 11445076 | SU11445076 | |
| Aesculap (old) | 3 metres | 11443043 | SU11443043 | |
| Aesculap (old) | 5 metres | 11445043 | SU11445043 | |

Bipolar cables coaxil pin

Coaxil pin for use with the following generators:
Erbe, Downs and Siemens

| Product | Length | Code reusable | Code single use | Ends |
|---------|--------|------------------|--------------------|------|
|---------|--------|------------------|--------------------|------|













COAXIL PIN

| Product | Length | Code reusable | Code single use | Ends |
|----------------|----------|------------------|--------------------|------|
| Euro | 3 metres | 11423075 | SU11423075 | |
| Euro | 5 metres | 11425075 | SU11425075 | |
| Storz | 3 metres | 11423041 | SU11423041 | |
| Storz | 5 metres | 11425041 | SU11425041 | |
| Wolf | 3 metres | 11423042 | SU11423042 | |
| Wolf | 5 metres | 11425042 | SU11425042 | |
| US | 3 metres | 11423076 | SU11423076 | |
| US | 5 metres | 11425076 | SU11425076 | |
| Aesculap (old) | 3 metres | 11423043 | SU11423043 | |
| Aesculap (old) | 5 metres | 11425043 | SU11425043 | |


Bipolar cables coaxial pin















Coaxial pin for use with the following generators:
Berchtold, Wolf and Martin

| Product | Length | Code reusable  | Code single use  | Ends |
|--------------------|----------|--|--|--|
| COAXIAL PIN | | | | |
| Euro | 3 metres | 11453075 | SU11453075 |  |
| Euro | 5 metres | 11455075 | SU11455075 |  |
| Storz | 3 metres | 11453041 | SU11453041 |  |
| Storz | 5 metres | 11455041 | SU11455041 |  |
| Wolf | 3 metres | 11453042 | SU11453042 |  |
| Wolf | 5 metres | 11455042 | SU11455042 |  |
| US | 3 metres | 11453076 | SU11453076 |  |
| US | 5 metres | 11455076 | SU11455076 |  |
| Aesculap (old) | 3 metres | 11453043 | SU11453043 |  |
| Aesculap (old) | 5 metres | 11455043 | SU11455043 |  |

Bipolar cables lemo connection



Lemo connection for use with the following generator:
Older Eshmann (lemo)

| Product | Length | Code reusable  | Code single use  | Ends |
|------------------------|----------|--|--|--|
| LEMO CONNECTION | | | | |
| Euro | 3 metres | 11433075 | SU11433075 |  |
| Euro | 5 metres | 11435075 | SU11435075 |  |
| Storz | 3 metres | 11433041 | SU11433041 |  |
| Storz | 5 metres | 11435041 | SU11435041 |  |
| Wolf | 3 metres | 11433042 | SU11433042 |  |
| Wolf | 5 metres | 11435042 | SU11435042 |  |
| US | 3 metres | 11433076 | SU11433076 |  |
| US | 5 metres | 11435076 | SU11435076 |  |
| Aesculap (old) | 3 metres | 11433043 | SU11433043 |  |
| Aesculap (old) | 5 metres | 11435043 | SU11435043 |  |

Bipolar cables fixed 4mm pins



Fixed 4mm pins for use with the following generators:
EMC and Dolley

| Product | Length | Code <small>reusable</small> | Code <small>single use</small> | Ends |
|---------|--------|---------------------------------|-----------------------------------|------|
|---------|--------|---------------------------------|-----------------------------------|------|

FIXED 4mm PINS

| | | | | |
|----------------|----------|----------|------------|--|
| Euro | 3 metres | 11483075 | SU11483075 | |
| Euro | 5 metres | 11485075 | SU11485075 | |
| Storz | 3 metres | 11483041 | SU11483041 | |
| Storz | 5 metres | 11485041 | SU11485041 | |
| Wolf | 3 metres | 11483042 | SU11483042 | |
| Wolf | 5 metres | 11485042 | SU11485042 | |
| US | 3 metres | 11483076 | SU11483076 | |
| US | 5 metres | 11485076 | SU11485076 | |
| Aesculap (old) | 3 metres | 11483043 | SU11483043 | |
| Aesculap (old) | 5 metres | 11485043 | SU11485043 | |

Bipolar cables fixed 4mm pins




Fixed 4mm pins for use with the following generators:
Valleylab and Lamidey

| Product | Length | Code <small>reusable</small> | Code <small>single use</small> | Ends |
|---------|--------|---------------------------------|-----------------------------------|------|
|---------|--------|---------------------------------|-----------------------------------|------|

FIXED 4mm PINS

| | | | | |
|----------------|----------|----------|------------|--|
| Euro | 3 metres | 11473075 | SU11473075 | |
| Euro | 5 metres | 11475075 | SU11475075 | |
| Storz | 3 metres | 11473041 | SU11473041 | |
| Storz | 5 metres | 11475041 | SU11475041 | |
| Wolf | 3 metres | 11473042 | SU11473042 | |
| Wolf | 5 metres | 11475042 | SU11475042 | |
| US | 3 metres | 11473076 | SU11473076 | |
| US | 5 metres | 11475076 | SU11475076 | |
| Aesculap (old) | 3 metres | 11473043 | SU11473043 | |
| Aesculap (old) | 5 metres | 11475043 | SU11475043 | |

Reusable bipolar forceps



- 1 Straight standard pattern
- 2 Straight cawthorne pattern
- 3 Straight adson pattern
- 4 Angled standard pattern
- 5 Bayonet standard pattern

| Product | Length | Tip size 2mm | Tip size 1mm | Tip size 0.5mm |
|---------|--------|-----------------|-----------------|-------------------|
|---------|--------|-----------------|-----------------|-------------------|

STRAIGHT FORCEPS STANDARD STYLE

| | | | | |
|------------------|-------|-----------|-----------|-----------|
| Standard pattern | 110mm | 7720110UK | 7710110UK | 7705110UK |
| Standard pattern | 130mm | 7720130UK | 7710130UK | 7705130UK |
| Standard pattern | 150mm | 7720150UK | 7710150UK | 7705150UK |
| Standard pattern | 180mm | 7720180UK | 7710180UK | 7705180UK |
| Standard pattern | 200mm | 7720200UK | 7710200UK | 7705200UK |

STRAIGHT FORCEPS CAWTHORNE STYLE

| | | | | |
|------------------|-------|-----------|-----------|-----------|
| Standard pattern | 110mm | 7720110UK | 7710110UK | 7705110UK |
| Standard pattern | 130mm | 7720130UK | 7710130UK | 7705130UK |

STRAIGHT FORCEPS ADSON STYLE

| | | | | |
|---------------|-------|-----------|-----------|-----------|
| Adson pattern | 130mm | 7920130UK | 7910130UK | 7905130UK |
| Adson pattern | 150mm | 7920150UK | 7910150UK | 7905150UK |

ANGLED FORCEPS STANDARD AND CAWTHORNE STYLE


| | | | | |
|------------------|-------|-----------|-----------|-----------|
| Standard pattern | 80mm | 7520080UK | 7510080UK | 7505080UK |
| Standard pattern | 110mm | 7520110UK | 7510110UK | 7505110UK |
| Standard pattern | 130mm | 7520130UK | 7510130UK | 7505130UK |
| Standard pattern | 150mm | 7520150UK | 7510150UK | 7505150UK |
| Standard pattern | 180mm | 7520180UK | 7510180UK | 7505180UK |
| Standard pattern | 200mm | 7520200UK | 7510200UK | 7505200UK |
| Standard pattern | 230mm | 7520230UK | 7510230UK | 7505230UK |

STRAIGHT FORCEPS BAYONET STYLE

| | | | | |
|-----------------|-------|-----------|-----------|-----------|
| Bayonet pattern | 150mm | 7820150UK | 7810150UK | 7805150UK |
| Bayonet pattern | 200mm | 7820200UK | 7810200UK | 7805200UK |

All codes shown are standard UK connection. US connection available please contact us for more details.

Single use bipolar and INSU-TIP forceps



Standard single use

| Product | Length | Tip size 2mm | Tip size 1mm | Tip size 0.5mm |
|---------|--------|-----------------|-----------------|-------------------|
|---------|--------|-----------------|-----------------|-------------------|

SINGLE USE STRAIGHT FORCEPS STANDARD STYLE

| | | | | |
|------------------------|-------|--------------|--------------|--------------|
| Single use standard | 110mm | SU772011044 | SU771011044 | SU770511044 |
| Single use standard | 130mm | SU772013044 | SU771013044 | SU770513044 |
| Single use standard | 150mm | SU772015044 | SU771015044 | SU770515044 |
| Single use standard | 180mm | SU772018044 | SU771018044 | SU770518044 |
| Single use standard | 200mm | SU772020044 | SU771020044 | SU770520044 |
| Single use guiding pin | 200mm | SU772020044P | SU771020044P | SU770520044P |
| INSU-TIP standard | 110mm | SU872011044 | SU871011044 | SU870511044 |
| INSU-TIP standard | 130mm | SU872013044 | SU871013044 | SU870513044 |
| INSU-TIP standard | 150mm | SU872015044 | SU871015044 | SU870515044 |
| INSU-TIP standard | 180mm | SU872018044 | SU871018044 | SU870518044 |
| INSU-TIP standard | 200mm | SU872020044 | SU871020044 | SU870520044 |

SINGLE USE FORCEPS CAWTHORNE STYLE

| | | | | |
|----------------------|-------|-------------|-------------|-------------|
| Single use cawthorne | 80mm | SU762008044 | SU761008044 | SU760508044 |
| Single use cawthorne | 110mm | SU762011044 | SU761011044 | SU760511044 |
| INSU-TIP cawthorne | 80mm | SU862008044 | SU861008044 | SU860508044 |
| INSU-TIP cawthorne | 110mm | SU862011044 | SU861011044 | SU860511044 |

SINGLE USE FORCEPS STANDARD ANGLED TIPS

| | | | | |
|---------------------|-------|-------------|-------------|-------------|
| Single use standard | 80mm | SU752008044 | SU751008044 | SU750508044 |
| Single use standard | 110mm | SU752011044 | SU751011044 | SU750511044 |
| Single use standard | 130mm | SU752013044 | SU751013044 | SU750513044 |
| Single use standard | 150mm | SU752015044 | SU751015044 | SU750515044 |
| Single use standard | 180mm | SU752018044 | SU751018044 | SU750518044 |
| Single use standard | 200mm | SU752020044 | SU751020044 | SU750520044 |
| INSU-TIP standard | 80mm | SU852008044 | SU851008044 | SU850508044 |
| INSU-TIP standard | 110mm | SU852011044 | SU851011044 | SU850511044 |
| INSU-TIP standard | 130mm | SU852013044 | SU851013044 | SU850513044 |
| INSU-TIP standard | 150mm | SU852015044 | SU851015044 | SU850515044 |
| INSU-TIP standard | 180mm | SU852018044 | SU851018044 | SU850518044 |
| INSU-TIP standard | 200mm | SU852020044 | SU851020044 | SU850520044 |

Single use bipolar and INSU-TIP forceps

| Product | Length | Tip size 2mm | Tip size 1mm | Tip size 0.5mm |
|---------|--------|-----------------|-----------------|-------------------|
|---------|--------|-----------------|-----------------|-------------------|

SINGLE USE FORCEPS BAYONET STYLE

| | | | | |
|--------------------|-------|-------------|-------------|-------------|
| Disposable bayonet | 150mm | SU782015044 | SU781015044 | SU780515044 |
| Disposable bayonet | 200mm | SU782020044 | SU781020044 | SU780520044 |
| INSU-TIP bayonet | 150mm | SU882015044 | SU881015044 | SU880515044 |
| INSU-TIP bayonet | 200mm | SU882020044 | SU881020044 | SU880520044 |

SINGLE USE FORCEPS ADSON STYLE

| | | | | |
|------------------|-------|-------------|-------------|-------------|
| Single use adson | 130mm | SU792013044 | SU791013044 | SU790513044 |
| Single use adson | 150mm | SU792015044 | SU791015044 | SU790515044 |
| INSU-TIP adson | 130mm | SU892013044 | SU891013044 | SU890513044 |
| INSU-TIP adson | 150mm | SU892015044 | SU891015044 | SU890515044 |

Single use and reusable monopolar forceps straight forceps



- 1 Straight standard pattern
- 2 Straight standard pattern with cable
- 3 Angled standard pattern
- 4 Bayonet standard pattern
- 5 Straight adson style
- 6 Straight McIndoe style


| Product | Length | Tip size 2mm | Tip size 1mm | Cable 2mm | Cable 1mm |
|---------|--------|-----------------|-----------------|--------------|--------------|
|---------|--------|-----------------|-----------------|--------------|--------------|

REUSABLE AND SINGLE USE FORCEPS VARIOUS STYLES

| | | | | | |
|---------------------|-------|-----------|-----------|-------------|-------------|
| Standard pattern | 150mm | 672015048 | 671015048 | SU672015040 | SU671015040 |
| Standard pattern | 180mm | 672018048 | 671018048 | SU672018040 | SU671018040 |
| Standard pattern | 200mm | 672020048 | 671020048 | SU672020040 | SU671020040 |
| Angled standard | 180mm | 652018048 | 651018048 | SU652018040 | SU651018040 |
| Bayonet straight | 150mm | 682015048 | 681015048 | SU682015040 | SU681015040 |
| Bayonet straight | 200mm | 682020048 | 681020048 | SU682020040 | SU681020040 |
| RemingtonH straight | 230mm | 632023048 | 631015048 | SU632015040 | SU631015040 |
| Adson straight | 130mm | 692013048 | 691013048 | SU692013040 | SU691013040 |
| McIndoe straight | 150mm | 642015048 | 641015048 | SU642015040 | SU641015040 |

Other sizes available upon request. Please contact us for more details

Monopolar scissors



| Product | Length | Straight <input type="checkbox"/> | Curved <input type="checkbox"/> |
|---------|--------|-----------------------------------|---------------------------------|
|---------|--------|-----------------------------------|---------------------------------|

SCISSOR STYLES REUSABLE

| | | | |
|----------------|-------|--------|--------|
| Turner warwick | 180mm | 601800 | 601801 |
| Turner warwick | 240mm | 602400 | 602401 |
| Wilson hey | 130mm | 611300 | 611301 |
| Wilson hey | 150mm | 611500 | 611501 |
| Wilson hey | 180mm | 611800 | 611801 |
| Wilson hey | 200mm | 612000 | 612001 |
| Wilson hey | 230mm | 612300 | 612301 |
| Riches | 130mm | 621300 | 621301 |
| Riches | 150mm | 621500 | 621501 |
| Riches | 180mm | 621800 | 621801 |
| Riches | 200mm | 622000 | 622001 |
| Riches | 230mm | 622300 | 622301 |

Single use available upon request. Please contact us for more details

1 Turner warwick 2 Wilson hey curved 3 Wilson hey straight 4 Riches

Electrodes laparoscopic



End options

| Product | Codes | Single use <input type="checkbox"/> |
|---------|-------|-------------------------------------|
|---------|-------|-------------------------------------|

MONOPOLAR LAPAROSCOPIC HANDLES ANGLED AND STRAIGHT STYLE

| | | |
|--|----------|------------|
| Monopolar laparoscopic handle angled | 3840201A | SU3840201A |
| Monopolar laparoscopic handle angled | 3840202A | SU3840202A |
| Monopolar laparoscopic handle angled | 3840203A | SU3840203A |
| Monopolar laparoscopic handle angled | 3840204A | SU3840204A |
| Monopolar laparoscopic handle straight | 3840201 | SU3840201 |
| Monopolar laparoscopic handle straight | 3840202 | SU3840202 |
| Monopolar laparoscopic handle straight | 3840203 | SU3840203 |
| Monopolar laparoscopic handle straight | 3840204 | SU3840204 |

| Product | Code | Box qty |
|---------|------|---------|
|---------|------|---------|

INSUFFLATORS

| | | |
|--|---------|----|
| Insufflator 15mm machine connection (WOLF) | 8870101 | 30 |
| Insufflator male luer connection (STORZ/OLYMPUS) | 8870900 | 30 |
| Insufflator female luer connection (GENERIC) | 8870901 | 30 |

Insufflators



| Product | Code | Box qty |
|---------|------|---------|
|---------|------|---------|

ENT INSTRUMENTS

| | | |
|----------------------------|--------|---|
| Abbey sub-mucosal needle | 885343 | 1 |
| Simpson pistol grip handle | 884256 | 1 |
| Killian nasal speculum | 884107 | 1 |

THUDICUMS

| | | | |
|-------------------------------------|-----------|----|---|
| Single use thudicum speculum size 1 | SU41110N1 | 96 | 2 |
| Single use thudicum speculum size 2 | SU41120N1 | 96 | 2 |
| Single use thudicum speculum size 3 | SU41130N1 | 96 | 2 |
| Single use thudicum speculum size 4 | SU41140N1 | 96 | 2 |
| Single use thudicum speculum size 5 | SU41150N1 | 96 | 2 |
| Single use thudicum speculum size 6 | SU41160N1 | 96 | 2 |
| Single use thudicum speculum size 7 | SU41170N1 | 96 | 2 |


DIATHERMY QUIVERS

| | | | |
|--|-----------|----|---|
| Quiver 65mm dia x 400m length telescopic | 1565400 | 1 | |
| Quiver 40mm dia x 235m length | 1540235 | 1 | |
| Quiver 65mm dia x 190m length | 1565190 | 1 | |
| Quiver 65mm dia x 235m length telescopic | 156235 | 1 | |
| Single use quiver 65mm dia x 190m length | SU1565190 | 96 | 2 |

Re-insulation
repair service
Insulation
thickness testing

Re-insulation
repair service
Re-insulation
respray service

Re-insulation
repair service
Insulation
pinhole testing



1 Abbey sub-mucosal 2 Simpson pistol grip handle 3 Killian nasal speculum 4 Thudicum 5 Diathermy telescopic quiver 65mm diameter x 400mm length 6 Diathermy quiver 7 Diathermy quiver

Thudicums also available as reusable. Please contact us for more details.

Loop and ball electrodes

Product

Codes

2.4mm



4mm



STANDARD FRAME

| | | |
|--|--------------|--------------|
| Standard frame 130mm 2.4mm shaft 08mm loop | SU3024130L08 | SU3040130L08 |
| Standard frame 130mm 2.4mm shaft 10mm loop | SU3024130L10 | SU3040130L10 |
| Standard frame 130mm 2.4mm shaft 12mm loop | SU3024130L12 | SU3040130L12 |
| Standard frame 130mm 2.4mm shaft 15mm loop | SU3024130L15 | SU3040130L15 |
| Standard frame 130mm 2.4mm shaft 18mm loop | SU3024130L18 | SU3040130L18 |
| Standard frame 130mm 2.4mm shaft 20mm loop | SU3024130L20 | SU3040130L20 |
| Standard frame 130mm 2.4mm shaft 22mm loop | SU3024130L22 | SU3040130L22 |
| Standard frame 130mm 2.4mm shaft 25mm loop | SU3024130L25 | SU3040130L25 |

Product

Codes

2.4mm



4mm



LOW PROFILE FRAME

| | | |
|---|---------------|---------------|
| Low profile 130mm 2.4mm shaft 08mm loop | SU3024130LP08 | SU3040130LP08 |
| Low profile 130mm 2.4mm shaft 10mm loop | SU3024130LP10 | SU3040130LP10 |
| Low profile 130mm 2.4mm shaft 12mm loop | SU3024130LP12 | SU3040130LP12 |
| Low profile 130mm 2.4mm shaft 15mm loop | SU3024130LP15 | SU3040130LP15 |
| Low profile 130mm 2.4mm shaft 18mm loop | SU3024130LP18 | SU3040130LP18 |
| Low profile 130mm 2.4mm shaft 20mm loop | SU3024130LP20 | SU3040130LP20 |
| Low profile 130mm 2.4mm shaft 22mm loop | SU3024130LP22 | SU3040130LP22 |
| Low profile 130mm 2.4mm shaft 25mm loop | SU3024130LP25 | SU3040130LP25 |

Product

Codes

2.4mm



4mm



Available in 3mm-8mm balls

BALL ELECTRODES

| | | |
|---|-------------|-------------|
| Ball electrode 130mm 2.4mm shaft 3mm ball | SU3024130B3 | SU3040130B3 |
| Ball electrode 130mm 2.4mm shaft 4mm ball | SU3024130B4 | SU3040130B4 |
| Ball electrode 130mm 2.4mm shaft 5mm ball | SU3024130B5 | SU3040130B5 |
| Ball electrode 130mm 2.4mm shaft 6mm ball | SU3024130B6 | SU3040130B6 |
| Ball electrode 130mm 2.4mm shaft 8mm loop | SU3024130B8 | SU3040130B8 |

Product

Codes

2.4mm



4mm



SAIL FRAME

| | | |
|--|---------------|---------------|
| Sail frame 130mm 2.4mm shaft 12mm loop | SU3024130SA12 | SU3040130SA12 |
| Sail frame 130mm 2.4mm shaft 15mm loop | SU3024130SA15 | SU3040130SA15 |
| Sail frame 130mm 2.4mm shaft 20mm loop | SU3024130SA20 | SU3040130SA20 |

Product

Codes

2.4mm



4mm



SQUARE FRAME

| | | |
|--|---------------|---------------|
| Square frame 130mm 2.4mm shaft 08mm loop | SU3024130SQ08 | SU3040130SQ08 |
| Square frame 130mm 2.4mm shaft 10mm loop | SU3024130SQ10 | SU3040130SQ10 |
| Square frame 130mm 2.4mm shaft 12mm loop | SU3024130SQ12 | SU3040130SQ12 |
| Square frame 130mm 2.4mm shaft 15mm loop | SU3024130SQ15 | SU3040130SQ15 |
| Square frame 130mm 2.4mm shaft 18mm loop | SU3024130SQ18 | SU3040130SQ18 |
| Square frame 130mm 2.4mm shaft 20mm loop | SU3024130SQ20 | SU3040130SQ20 |
| Square frame 130mm 2.4mm shaft 22mm loop | SU3024130SQ22 | SU3040130SQ22 |
| Square frame 130mm 2.4mm shaft 25mm loop | SU3024130SQ25 | SU3040130SQ25 |

Product

Codes

2.4mm



4mm



T-BAR FRAME

| | | |
|---|--------------|--------------|
| T-bar frame 130mm 2.4mm shaft 12mm loop | SU3024130T12 | SU3040130T12 |
| T-bar frame 130mm 2.4mm shaft 15mm loop | SU3024130T15 | SU3040130T15 |
| T-bar frame 130mm 2.4mm shaft 20mm loop | SU3024130T20 | SU3040130T20 |

Product

Codes

2.4mm



4mm



TRAFFORD FRAME

| | | |
|--|-----------------|-----------------|
| Trafford frame 130mm 2.4mm shaft 10mm loop | SU3024130TR10 | SU3040130TR10 |
| Trafford frame 130mm 2.4mm shaft 15mm x 8mm | SU3024130TR1508 | SU3040130TR1508 |
| Trafford frame 130mm 2.4mm shaft 18mm x 12mm | SU3024130TR1812 | SU3040130TR1812 |

SMOKE EVACUATION KITS

An all in one filter containing a standard hepa with Carbon filtration.

Meditechs' all in one single use hepa carbon filter will reduce and eradicate residual smells and contaminants. The smoke is captured by the carbon hepa filter before it enters the smoke evacuator ensuring safety and continuity by single use application. Pre-assembled systems include:

- Filter elephant tubing connector.
- Cusco connector tubing.
- Connectors allow fitting to Erbe, ValleyLab, Conmed and other well known smoke evacuator machines.

Smoke evacuation pencil



Smoke pencil with tubing set SU37407M smoke pencil with tubing set and filter SU374075



Speculums




| Product | Size | Codes (1) (insulated) | Codes (1) (non-insulated) | Codes (1) (insulated) | Codes (1) (non-insulated) | Codes (1) (Plastic) |
|---------|------|-----------------------|---------------------------|-----------------------|---------------------------|---------------------|
|---------|------|-----------------------|---------------------------|-----------------------|---------------------------|---------------------|

GYNAECOLOGY INSTRUMENTS SPECULUM DESIGNS


| | | | | | | |
|-----------|-------------|--------|--------|-------------|--------------|-------------|
| Cusco | SMALL | 32010* | 32010N | SU32010 (2) | SU32010N (2) | SU61101 (2) |
| Cusco | MEDIUM | 32011* | 32011N | SU32011 (2) | SU32011N (2) | SU61102 (2) |
| Cusco | LARGE | 32012* | 32012N | SU32012 (2) | SU32012N (2) | SU61103 (2) |
| Cusco | MEDIUM LONG | 32013* | 32013N | SU32013 (2) | SU32013N (2) | |
| Cusco | EXTRA LARGE | 32016* | 32016N | SU32016 (2) | SU32016N (2) | |
| Cusco | MIXED | | | | SU611003 (2) | |
| Sims | SMALL | 33010 | 33010N | SU33010 (2) | SU33010N (2) | |
| Sims | MEDIUM | 33011 | 33011N | SU33011 (2) | SU33011N (2) | |
| Sims | LARGE | 33012 | 33012N | SU33012 (2) | SU33012N (2) | |
| Graves | SMALL | 34010 | 34010N | SU34010 (2) | SU34010N (2) | |
| Graves | MEDIUM | 34011 | 34011N | SU34011 (2) | SU34011N (2) | |
| Graves | LARGE | 34012 | 34012N | SU34012 (2) | SU34012N (2) | |
| Langden | SMALL | 35010 | 35010 | SU35010 (2) | SU35010N (2) | |
| Langden | MEDIUM | 35011 | 35011 | SU35011 (2) | SU35011N (2) | |
| Langden | LARGE | 35012 | 35012 | SU35012 (2) | SU35012N (2) | |
| Winterton | EXTRA LARGE | 32016 | 35016 | SU35016 (2) | SU35016N (2) | |

Disposable kits for smoke evacuation available with or without hepa filter. Portable smoke evacuator also available.

Biopsy instruments



The Tishler Morgan is made of a high quality stainless steel with a precision biopsy head and is a very durable product.



- 1 Tishler Morgan with teeth close up
- 2 Eppendorfer with teeth close up
- 3 Rotational biopsy punch

| Product | Code | Single use code |
|---------|------|-----------------|
|---------|------|-----------------|

TISHLER MORGAN BIOPSY INSTRUMENT

| | | |
|----------------------------------|---------|-----------|
| Tishler morgan biopsy instrument | 5720021 | SU5720021 |
|----------------------------------|---------|-----------|

| Product | Code | Single use code |
|---------|------|-----------------|
|---------|------|-----------------|

EPPENDORFER BIOPSY INSTRUMENT

| | | |
|--|---------|-----------|
| Eppendorfer biopsy instrument 200mm length | 5720022 | SU5720022 |
|--|---------|-----------|

| | | |
|--|---------|-----------|
| Eppendorfer biopsy instrument 220mm length | 5722022 | SU5722022 |
|--|---------|-----------|

| Product | Code | Single use code |
|---------|------|-----------------|
|---------|------|-----------------|

EPPENDORFER BIOPSY INSTRUMENT



| | | |
|---|---------|-----------|
| Tishler rotational biopsy punch 200mm shaft | 5721120 | SU5721120 |
|---|---------|-----------|

| | | |
|--|---------|-----------|
| Morgan rotational biopsy punch 220mm shaft | 5721122 | SU5721122 |
|--|---------|-----------|


| | | |
|--|---------|-----------|
| Mini rotational biopsy punch 200mm shaft | 5721220 | SU5721220 |
|--|---------|-----------|

| | | |
|---|---------|-----------|
| Tishler rotational biopsy punch 220mm shaft | 5721222 | SU5721222 |
|---|---------|-----------|

Single use consumables



 AVAILABLE
AS REUSABLE
IF REQUIRED

- 1 Endocervical brush
- 2 Amniotic membrane hook
- 3 Cervical smear wood spatula
- 4 Endometrial sampler uterine mucosal
- 5 Cell brush
- 6 Single use double ended speculum (dilatator) 5/6 220mm
- 7 Sponge forceps plastic
- 8 Vulsellum forceps metal

 NEW
NEEDLE SAFE
SYSTEM
AVAILABLE


| Product | Code | Box qty |
|---|---------------|---------|
| CONSUMABLE PRODUCTS | | |
| Endocervical brush | 720105 | 100 |
| Single use plastic amniotic membrane hook | 720193 | 100 |
| Cervical smear wood spatula | 720194 | 100 |
| Endometrial sampler uterine mucosal | 7201935 | 100 |
| Cell brush | 720196 | 100 |
| Single use sponge forceps metal | SU57220031 | 10 |
| Single use sponge forceps plastic | SU57220031P | 10 |
| Single use vulsellum forceps metal | SU57229900 | 20 |
| Single use vulsellum forceps plastic | SU57229900P | 20 |
| Single use dental syringe | SU57229920 | 20 |
| Single use dental syringe long | SU57229920L | 20 |
| Single use dental syringe extra long side load | SU57229920XL | 20 |
| Single use tenaculum | SU572TEN | 20 |
| Single use polypus/ovum forceps 250mm | SUB811166 | 10 |
| Single use kogans 9 1/2 calibrated | SUB816KG | 10 |
| Single use kurihara 250mm endocervical speculae | SUB8KH250 | 10 |
| SU double ended speculum (dilatators) 1/2 220mm | SUB8DESP22015 | 20 |
| SU double ended speculum (dilatators) 3/4 220mm | SUB8DESP22017 | 20 |
| SU double ended speculum (dilatators) 5/6 220mm | SUB8DESP22021 | 20 |

Theatre and clinic consumables




- 1 Iris scissors
- 2 Metzenbaum scissors
- 3 Dressing scissors
- 4 Mayo scissors
- 5 Stitch scissors
- 6 Spencer wells forceps
- 7 Supercut scissors

| Product | Single use code | Box qty |
|--|-----------------|---------|
| THEATRE AND CLINIC CONSUMABLE PRODUCTS | | |
| Single use iris scissors 11.5cm | SU88AB23211K | 20 |
| Single use iris scissors 11.5cm curved | SU88AB23211C | 20 |
| Single use metzenbaum scissors 23cm | SU88B11523K | 20 |
| Single use metzenbaum scissors 23cm curved | SU88B11523CK | 20 |
| Single use dressing scissors 13cm blunt/blunt | SU88B8BL130K | 20 |
| Single use dressing scissors 13cm blunt/sharp | SU88PBS130 | 20 |
| Single use dressing scissors sharp/sharp | SU88SSD | 20 |
| Single use mayo scissors 170mm | SU88SMA170K | 20 |
| Single use mayo scissors 170mm curved | SU88SMA170CK | 20 |
| Single use stitch scissors 5.5 | SU88SSTK | 20 |
| Single use stitch scissors 4.4 Sharp | SU88ST115K | 20 |
| Single use scissors angled plastic handle supercut | SU88STCK | 20 |
| Single use lister bandage scissors | SU88BDLK | 20 |
| Single use umbilical scissors 100mm | SU88U100 | 20 |
| Single use spencer wells forcep 125cm | SU88ASW125 | 20 |
| Single use spencer wells forcep 140mm | SU88ASW140 | 20 |
| Single use spencer wells forcep 180mm | SU88ASW180 | 20 |
| Single use spencer wells forcep 190mm | SU88ASW190 | 20 |
| Single use spencer wells forcep 230mm | SU88ASW230 | 20 |

Retractable safety scalpels

1



| Blade size | Code | Blade | Code |
|------------|------|-------|------|
|------------|------|-------|------|

UNIVERSAL RETRACTABLE SAFETY SCALPELS

| | | | |
|-----------|-----------|-----------|-----------|
| Blade 9 | SUBBSS09 | Blade 17 | SU88SS17 |
| Blade 10 | SUBBSS10 | Blade 18 | SU88SS18 |
| Blade 10A | SUBBSS10A | Blade 19 | SU88SS19 |
| Blade 11 | SUBBSS11 | Blade 20 | SU88SS20 |
| Blade 11K | SUBBSS11K | Blade 21 | SU88SS21 |
| Blade 13 | SUBBSS13 | Blade 22 | SU88SS22 |
| Blade 14 | SUBBSS14 | Blade 22A | SU88SS22A |
| Blade 15 | SUBBSS15 | Blade 23 | SU88SS23 |
| Blade 15B | SUBBSS15B | Blade 24 | SU88SS24 |
| Blade 15C | SUBBSS15C | Blade 25 | SU88SS25 |
| Blade 16 | SUBBSS16 | | |

Stainless steel handles

| Blade size | Code | Product |
|------------|------|---------|
|------------|------|---------|

METAL HANDLES

| | | |
|-----------|--------------|--|
| Blade 9 | PSSH3/3L/7/9 | |
| Blade 10 | PSSH3/3L/7/9 | |
| Blade 10A | PSSH3/3L/7/9 | |
| Blade 11 | PSSH3/3L/7/9 | |
| Blade 11K | PSSH3/3L/7/9 | |
| Blade 13 | PSSH3/3L/7/9 | |
| Blade 14 | PSSH3/3L/7/9 | |
| Blade 15 | PSSH3/3L/7/9 | |
| Blade 15B | PSSH3/3L/7/9 | |
| Blade 15C | PSSH3/3L/7/9 | |
| Blade 16 | PSSH3/3L/7/9 | |
| Blade 17 | PSSH3/3L/7/9 | |
| Blade 18 | PSSH4/4L/6 | |
| Blade 19 | PSSH4/4L/6 | |
| Blade 20 | PSSH4/4L/6 | |
| Blade 21 | PSSH4/4L/6 | |
| Blade 22 | PSSH4/4L/6 | |
| Blade 22A | PSSH4/4L/6 | |
| Blade 23 | PSSH4/4L/6 | |
| Blade 24 | PSSH4/4L/6 | |
| Blade 25 | PSSH4/4L/6 | |

Surgical blades

| Blade size | Code | Product |
|------------|------|---------|
|------------|------|---------|

SURGICAL BLADES

| | | |
|-----------|--------|--|
| Blade 9 | PSS09 | |
| Blade 10 | PSS10 | |
| Blade 10A | PSS10A | |
| Blade 11 | PSS11 | |
| Blade 11K | PSS11K | |
| Blade 12 | PSS12 | |
| Blade 13 | PSS13 | |
| Blade 14 | PSS14 | |
| Blade 15 | PSS15 | |
| Blade 15B | PSS15B | |
| Blade 15C | PSS15C | |
| Blade 16 | PSS16 | |
| Blade 17 | PSS17 | |
| Blade 18 | PSS18 | |
| Blade 19 | PSS19 | |
| Blade 20 | PSS20 | |
| Blade 21 | PSS21 | |
| Blade 22 | PSS22 | |
| Blade 22A | PSS22A | |
| Blade 23 | PSS23 | |
| Blade 24 | PSS24 | |
| Blade 25 | PSS25 | |

Handles shown from blades 9 – blade 17 are PSSH7. Handles shown from blades 18 – blade 25 are PSSH4L

Podiatry products



2

| Product | Single use code | Box qty |
|---------|-----------------|---------|
|---------|-----------------|---------|

PODIATRY PRODUCTS

| | | |
|---|------------|----|
| Single use scalpel D15/E11 size 10 (also 11 and 15) | SU88PS10* | 20 |
| Single use diamond deb file | SU88PDDF | 20 |
| Single use emery board | SU88PEB | 20 |
| Single use nippers 14cm | SU88PN140C | 20 |
| Single blacks file | SU88PBF | 20 |
| Single use dressing scissors sharp/sharp 15cm | SU88SSD | 20 |
| Single use dressing scissors 15cm blunt/sharp | SU88PBS150 | 20 |
| Single use nippers 15cm | SU88PN150 | 20 |
| Dressing scissors | SU88PBS130 | 20 |

*Change end digits to your size specification



3



4

- 1 Safety scalpel
- 2 Blacks file medium
- 3 Nail nipper 14cm
- 4 Emery board

| Product | Single use code |
|---------|-----------------|
|---------|-----------------|

PODIATRY PACKS

PODIATRY PACK 1 | BASIC FOOTCARE SET:

Includes: nail nipper 14cm, curved blacks file medium & emery board

SU88PED01

PODIATRY PACK 2 | NAIL SURGERY SET:

Includes: scalpel handle 13.5cm, nippers elevator 11.5cm, dressing scissors sharp/blunt 14cm, blacks file-fine.

SU88PED02

Skin biopsy punch and dermal curettes



FEATURES

- Ribbed handle for effective grip and control.
- Seamless stainless steel sharp cutting edge.
- Size embossed on each punch.
- Wide range of sizes.
- Sterile and ready for use.
- Disposable.

Product

Single use code

Box qty

SKIN BIOPSY PUNCH

| | | |
|------------------------------------|-------------|---|
| Single use skin biopsy punch 1mm | SUB8SBP1SS | ☒ |
| Single use skin biopsy punch 1.5mm | SUB8SBP15SS | ☒ |
| Single use skin biopsy punch 2mm | SUB8SBP2SS | ☒ |
| Single use skin biopsy punch 2.5mm | SUB8SBP25SS | ☒ |
| Single use skin biopsy punch 3mm | SUB8SBP3SS | ☒ |
| Single use skin biopsy punch 3.5mm | SUB8SBP35SS | ☒ |
| Single use skin biopsy punch 4mm | SUB8SBP4SS | ☒ |
| Single use skin biopsy punch 5mm | SUB8SBP5SS | ☒ |
| Single use skin biopsy punch 6mm | SUB8SBP6SS | ☒ |
| Single use skin biopsy punch 8mm | SUB8SBP8SS | ☒ |



Product

Single use code

Box qty

DERMAL CURETTES

| | | |
|-------------------------------|----------|---|
| Single use dermal curette 2mm | SUB8DC20 | ☒ |
| Single use dermal curette 3mm | SUB8DC30 | ☒ |
| Single use dermal curette 4mm | SUB8DC40 | ☒ |
| Single use dermal curette 5mm | SUB8DC50 | ☒ |
| Single use dermal curette 7mm | SUB8DC70 | ☒ |

FEATURES

- Seamless stainless steel sharp cutting edge.
- Easy dynamic gripping handle.
- Sterile and ready for use.
- Disposable.

Thudicums and fraizer suction products




| Product | Single use code | Box qty |
|---------|-----------------|---------|
|---------|-----------------|---------|

THUDICUMS

| | | |
|-------------------------------------|-----------|----|
| Single use thudicum speculum size 1 | SU41110N1 | 20 |
| Single use thudicum speculum size 2 | SU41120N1 | 20 |
| Single use thudicum speculum size 3 | SU41130N1 | 20 |
| Single use thudicum speculum size 4 | SU41140N1 | 20 |
| Single use thudicum speculum size 5 | SU41150N1 | 20 |
| Single use thudicum speculum size 6 | SU41160N1 | 20 |
| Single use thudicum speculum size 7 | SU41170N1 | 20 |

Thudicums also available as reusable. Please contact us for more details.



| Product | Single use code | Box qty |
|---------|-----------------|---------|
|---------|-----------------|---------|

SUCTION PRODUCTS

| | | |
|--|-----------|----|
| Poole suction handle | SU8832008 | 20 |
| Fraizer suction tip with malleable shaft & control vent 8FR | SU8832026 | 20 |
| Fraizer suction tip with malleable shaft & control vent 10FR | SU8832027 | 20 |
| Fraizer suction tip with malleable shaft & control vent 12FR | SU8832028 | 20 |



meditech Supporting every breath
systems limited

Meditech Systems Limited.
Shrublands Estate, Sherstock,
Shaftesbury, Dorset, SP7 9PT.

tel: +44 (0)1747 821546
web: meditechsystems.co.uk
email: sales@meditechsystems.co.uk