

Anaesthetic breathing systems


Introducing the Meditech range of anaesthetic breathing systems, also known as circuits, both standard and Antimicrobial/SilverSentinel.

The Meditech range of circuits are available in different lengths and configurations. Our range includes MEDIFLEX, which is a standard flexible tubing system, X-TENDIBLE, which is a length adjustable system and Smoothbore tubing systems. We also offer the MEDILIMB coaxial and TWIN TUBE single limb options where required. Mapleson configurations are also offered using an APL valve system adult and paediatric.


SILVERSENTINEL Antimicrobial systems

SilverSentinel systems are our range of Antimicrobial breathing circuits. They contain a silver ion additive within the tubing that rapidly disrupts the normal enzyme activities of certain bacteria resulting in rapidly combating Hospital Acquired Infections (HAI) such as MRSA and other micro-organisms. All SilverSentinel systems are identified with a distinct purple tinted tubing. SilverSentinel is available on all circuits except smoothbore.


Smoothbore heated wire circuit 1.8m.



Deluxe bain breathing circuits



Description	Code	Box qty
SMOOTHBORE CIRCUIT		
Adult Smoothbore Circuit 1.8m	22SB1822FF	10
SMOOTHBORE HEATED WIRE CIRCUIT		
Single Heated Wire Circuit 1.8m with 0.5m Limb	22SBHW1822MMU	10
X-TENDIBLE CIRCUIT		
Adult 1.8m Extendible Circuit	22XP1822FF	20
DELUXE BAIN CIRCUIT		
Deluxe Bain circuit 1.8m with APL valve (60cm) & 2L bag	22B1822APL	10




1. Smoothbore tubing circuit.
2. X-TENDIBLE tubing circuit.

A great space saver


QUALITY | VALUE | INTEGRITY

Anaesthesia breathing systems for GE Healthcare® Carestation®


Breathing systems for the Aisys®, Avance®, and Aespire® Carestation®



Breathing Systems with spirometry kit.



Breathing Systems for the ADU Carestation®



Customisation of Carestation® breathing systems available upon request. Please ask one of our representatives or visit us at www.meditechsystems.co.uk


*Aisys, Avance, Aespire, Carestation and ADU are all trademarks of GE Company.

TIMETAG™


A time management device which will reduce your manual recordings and record daily progress. Assists with change frequency, traceability and expiry.

Meditech Systems 7 day circuits now come with a unique TIMETAG Simply activate the tag with a finger press button. The red indicator highlights the daily progress. The unique reference number on the timer is linked to the product batch number in use. Other time control periods available ie: 48hrs, 24hrs please ask our representative for details.

- Visual time reference on product.
- Reduce manual date recording.



A 7 Day TIMETAG will be on your circuit for easy referencing see diagram opposite.




Active

1 day progress

7 day progress

Other time control periods are available upon request.



QUALITY | VALUE | INTEGRITY

Meditech Systems Ltd, Unit 3 Shrublands Estate, Sherstock, Shaftesbury, Dorset, SP7 9PT. UK
 tel: +44 (0)1747 821546 web: meditechsystems.co.uk email: sales@meditechsystems.co.uk